


PROTECTION AT THE SOURCE


SHIELDS FOR EVERY PROCEDURE




The ScatPad is a sterile scatter-radiation shield placed on top of the patient drape — a protective zone that reduces harmful exposure. Radiation dose is cumulative, and multiple low-level procedures each day add up for physicians, nurses & techs.

- ★ Reduces exposure time to physician & staff
- ★ Lightweight — conforms to patient anatomy
- ★ Sterile, single-use — no special disposal
- ★ Independently tested attenuation
- ★ Available in 75%, 90% & 95% at 90 kVp


COMMON SHIELD STYLES




MULTIPURPOSE
SMA 102



FEMORAL ANGIOGRAPHY
SAFA 301 / 302



EP · DEVICE IMPLANT
SEPA 401 / 402



EP · CENTER SCOOP
SEPA 403



PERIPHERAL
SPA 501 / 502

SCATPAD PART NUMBERS & SIZES

PART	DESCRIPTION	PROT	SIZE
SMA102	Multipurpose interventional	90%	12.5×16.5
SMA123	Multipurpose, large	95%	14.5×16.5
SAFA301	Femoral angio / heart cath	75%	14.5×16.5
SAFA302	Femoral angio / heart cath	90%	14.5×16.5
SEPA401	EP / cath device implant	75%	12.5×16.5
SEPA402-DT	EP / cath device implant	90%	12.5×16.5
SEPA403-CTRSCP	EP / cath, center scoop	95%	12.5×16.5

PART	DESCRIPTION	PROT	SIZE
SEPA412	EP / cath, large rectangular cutout	90%	12.5×16.5
SEPA453-X	Cath/EP BPA, XL w/ scoop	95%	16×28
SPA501	Peripheral	75%	11×34
SPA502	Peripheral	90%	11×34
SJAA601-X	Radial / jugular / TIPS access	75%	12.5×16.5
SDFA701	Dialysis access, de-clotting	75%	32×50
SDFA702	Dialysis access, de-clotting	90%	32×50



1 (800) 227-3633

cardiomedicalproducts.com · orders@cardiomedicalproducts.com
385 Franklin Ave, Suite L · Rockaway, NJ 07866

