



# SCATPAD™

## DISPOSABLE SCATTER RADIATION PROTECTION



- PROTECT YOUR SKULL AND BRAIN WITH THE HEADLINER
- LEAD FREE
- DISPOSABLE
- EXTRA BARRIER MATERIAL IN BOTH LENGTH AND HEIGHT LIGHTWEIGHT
- THREE PROTECTION LEVELS
- VENTED TOP

Headliner Cap		
Model	Description	Case Qty
CM-104	Scatter Radiation Protective Cap Attenuation level 75% .125 LE	20
CM-105	Scatter Radiation Protective Cap Attenuation level 90% .25 LE	20
CM-106	Scatter Radiation Protective Cap Attenuation level 95% .375 LE	20

ScatPad Part Numbers and Description			
Model	Description / Attenuation Levels / Case of 10	Dimensions	Procedure
SMA 101	Multipurpose Shield, Absorbent, No Fenestration 75%	12.5"X16.5"	Interventional
SMA 102	Multipurpose Shield, Absorbent, No Fenestration 90%	12.5"X16.5"	Interventional
SMA 123	Multipurpose Shield, Absorbent, No Fenestration 95%	14.5X16.5"	Interventional
SMFA 201	Multipurpose Shield, Absorbent, With Fenestration 75%	12.5"X16.5"	Biopsy
SMFA 202	Multipurpose Shield, Absorbent, With Fenestration 90%	12.5"X16.5"	Biopsy
SAFA 301	Angiography Shield, Absorbent, With Fenestration 75%	14.5"X16.5"	Femoral
SAFA 302	Angiography Shield, Absorbent, With Fenestration 90%	14.5"X16.5"	Femoral
SEPA 401-TT	Electrophysiology Shield, Absorbent, With Fenestration 75%	12.5"X16.5"	Left Sub, Cath/EP
SEPA 402-A	Electrophysiology Shield, Absorbent, With Fenestration 90%	12.5"X16.5"	Left Sub, Cath/EP
SEPA-402-DT	Electrophysiology Shield, Absorbent, With Fenestration 90%	12.5"X16.5"	Left/Right Cath/EP
SEPA 403-CTSCP	Electrophysiology Shield, Absorbent, With Rectangular Fenestration 95%	12.5"x16.5"	Left Sub, Cath/EP
SEPA 412	Electrophysiology Shield, Absorbent, With Rectangular Fenestration 90%	12"X20"	EP/ Ablation
SEPA 422	Electrophysiology Shield, Absorbent, With Scoop Fenestration, Extended 90%	14.5"x16.5"	Left Sub, Cath/EP
SEPA 453-X	EP/CATH Specialty X-Large Shield, Absorbent, With Fenestration 95%	16"X28"	EP/CATH
SPA 501	Peripheral Shield, Absorbent, Without Fenestration 75%	11"X34"	Peripheral
SPA 502	Peripheral Shield, Absorbent, Without Fenestration 90%	11"X34"	Peripheral
SDFA 701	Dialysis Access Shield, Absorbent 75%	32"x50"	Dialysis
SDFA 702	Dialysis Access Shield, Absorbent 90%	32"x50"	Dialysis



TO ORDER OR START YOUR TRIAL CONTACT US AT  
(800) 227-3633 OR INFO@CARDIOMEDICALPRODUCTS.COM