

Cardio Clear®

DEFIBRILLATION ELECTRODES

For defibrillation, synchronized cardioversion, external pacing & ECG monitoring



THE CARDIO CLEAR ADVANTAGE FEATURES & BENEFITS

- ★ Smaller pad with larger conductive surface area
- ★ Flexible material easily conforms with patient's anatomy
- ★ Two extended cord lengths — 60" and 120"
- ★ Available in both radiotransparent & radiotranslucent
- ★ Available in both adult & pediatric models
- ★ All major connector styles
- ★ Leads-Out / pre-connect models available
- ★ Latex-free
- ★ Extended shelf life
36 months



DISTRIBUTED IN BOXES OF 10

Ten electrodes per box



NOT FOR USE WITH AED DEFIBRILLATORS

CARDIO CLEAR · ORDERING GUIDE

CONNECTOR PART NUMBERS



ZOLL

Transparent · Translucent

Z100	Transparent
TZ200	Translucent
LTZ200	Translucent, Leads Out
LOCZ500	Transparent, Leads Out
PZ300	Pediatric
LZ100-LONG-LI	Transparent · 10 ft leads NEW
TZ200-LONG	Translucent · 10 ft leads NEW



PHYSIO-CONTROL

Medtronic / LIFEPAK

O101	Transparent
T0201	Translucent
LT0201	Translucent, Leads Out
LOC0501	Transparent, Leads Out
P0301	Pediatric
O101-LONG	Transparent · 10 ft leads NEW



PHILIPS — PLUG

Plug-style

P103	Transparent
TP203	Translucent
LTP203	Translucent, Leads Out
LOC503	Transparent, Leads Out
PP303	Pediatric



PHILIPS — BARREL

Barrel-style

H104	Transparent
TH204	Translucent
LTH204	Translucent, Leads Out



ANDERSON

Anderson-style

A102	Transparent
TA202	Translucent
LTA202	Translucent, Leads Out

OPTIONS & FIT — HOW TO READ THE CODES

TRANSPARENT

Radiotransparent — clear under fluoroscopy

LEADS OUT

Pre-connect; wires exit for fast set-up

LONG LEAD

10 ft wires for extended reach

TRANSLUCENT

Radiotranslucent — cost-effective

PEDIATRIC

Sized for pediatric patients

CONNECTORS

Zoll, Physio-Control (Medtronic), Philips & Anderson

ORDER NOW

1 (800) 227-3633

cardiomedicalproducts.com

orders@cardiomedicalproducts.com



CARDIO MEDICAL PRODUCTS, INC.

40 YEARS OF EXCELLENCE

385 Franklin Ave, Suite L, Rockaway, NJ 07866 · Made in the U.S.A.